



SPECIFICATION

Astragaloside II

Product Code:	OA09407
Synonyms:	(3b,6a,16b, 20R,24S)-3-[(2-O-Acetyl-b-D-xylopyranosyl)oxy]-20,24-Epoxy-16,25-dihydroxy-9,19-cyclolanostan-6-yl b-D-glucopyranoside Astrasieversianin VIII
CAS Number:	84676-89-1
Chemical Formula:	$C_{43}H_{70}O_{15}$
Molecular Weight:	827.01

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 98%
Water content (Karl Fischer):	max 3%