



SPECIFICATION

Raloxifene 4'-D-glucuronide

Product Code:	MR07277
Synonyms:	4-[6-Hydroxy-3-[4-2-[2-(1-piperidinyloxy)benzoyl]benzothien-2-yl]phenyl]-D-glucopyranosiduronic acid
CAS Number:	182507-22-8
Chemical Formula:	$C_{34}H_{35}NO_{10}S$
Molecular Weight:	649.71

TECHNICAL SPECIFICATION

Appearance:	Yellow crystalline powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure