



## SPECIFICATION

### 2-Nitrophenyl 2-acetamido-2-deoxy- $\alpha$ -D-glucopyranoside

Product Code:	EN03206
Synonyms:	ONP- $\alpha$ -D-GlcNAc 2-Nitrophenyl N-acetyl- $\alpha$ -D-glucosamine
CAS Number:	10139-01-2
Chemical Formula:	$C_{14}H_{18}N_2O_8$
Molecular Weight:	342.3

## TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity:	min 97%
Identity ( $^1H$ NMR):	Conforms to structure