



SPECIFICATION

D-Mandelonitrile-b-gentiobioside

Product Code:	OM05087
Synonyms:	D-Mandelonitrile-b-D-glucoside-6-O-b-D-glucoside D-Amygdalin α-[(6-O-b-D-Glucopyranosyl-b-D-glucopyranosyl)oxy]- (R)- Benzeneacetonitrile
CAS Number:	29883-15-6
Chemical Formula:	$C_{20}H_{27}NO_{11}$
Molecular Weight:	457.43

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 98%
Loss on drying:	max 0.5%
Heavy metals:	max 5ppm