



## SPECIFICATION

### 6-Chloro-3-indolyl b-D-cellobioside

Product Code:	EC03250
CAS Number:	425427-87-8
Chemical Formula:	$C_{20}H_{26}ClNO_{11}$
Molecular Weight:	491.87

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 98%
Specific optical rotation $[\alpha]_D^{20}$	-40 to -36° $^{20}$ (c=1, H <sub>2</sub> O):
Identity (IR):	Conforms to structure
Solubility (1% in water):	Clear, colourless solution
Water content (Karl Fischer):	max 5%