



## SPECIFICATION

### 1,3,4,6-Tetra-O-acetyl-2-O-trifluoromethanesulfonyl- $\beta$ -D-mannopyranose

Product Code:	MT06933
Synonyms:	1,3,4,6-Tetraacetate 2-(trifluoromethane-sulphonate)- $\beta$ -D-mannopyranose Mannose triflate TATM
CAS Number:	92051-23-5
Chemical Formula:	$C_{15}H_{19}F_3O_{12}S$
Molecular Weight:	480.37

## TECHNICAL SPECIFICATION

Appearance:	White to colourless crystalline powder
Purity (HPLC):	min 98%
Melting point:	118 - 122°C
Identity ( $^1H$ NMR):	Conforms to structure