



## SPECIFICATION

### Teniposide

Product Code:	FT09381
Synonyms:	4'-Dimethylepipodophyllotoxin-(4,6-O-2-thenylidene)-b-D-glucopyranoside
CAS Number:	29767-20-2
Chemical Formula:	$C_{32}H_{32}O_{13}S$
Molecular Weight:	656.65

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min. 98%
Specific optical rotation:	106 to 118°
Identity ( $^1H$ NMR):	Conforms to structure
Water content:	max.1.0%