



## SPECIFICATION

### Febuxostat

Product Code:	FF10814
Synonyms:	2-[3-Cyano-4-(2-methylpropoxy)phenyl]-4-methyl-5-thiazolecarboxylic acid
CAS Number:	144060-53-7
Chemical Formula:	$C_{16}H_{16}N_2O_3S$
Molecular Weight:	316.38

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 99%
Loss on drying:	max 1%
Residue on ignition:	max 0.1%
Heavy metals:	max 10ppm