



SPECIFICATION

L-Rhamnose monohydrate

Product Code:	MR00035
Synonyms:	6-Deoxy-L-mannose
CAS Number:	6155-35-7
Chemical Formula:	$C_6H_{12}O_5 \cdot H_2O$
Molecular Weight:	182.17

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 99%
Residue on Ignition:	max 0.2%
Heavy Metals:	max 5ppm
Calcium:	max 10ppm
Water content (Karl Fischer):	10%