



SPECIFICATION

Methyl 6-O-(N-heptylcarbamoyl)- α -D-glucopyranoside

Product Code:	DM06833
Synonyms:	6-O-(N-Heptylcarbamoyl)-methyl- α -D-glucopyranoside HECAMEG
CAS Number:	115457-83-5
Chemical Formula:	$C_{15}H_{29}NO_7$
Molecular Weight:	335.39

TECHNICAL SPECIFICATION

Appearance:	White crystalline powder
Purity (HPLC):	min 98%
Purity (TLC):	One spot
Solubility (10 % in water):	Clear, colourless liquid
Specific optical rotation $[\alpha]_D$	Reported 20 (c=1, H ₂ O):