



SPECIFICATION

Benzamidine HCl hydrate

Product Code:	FB05376
CAS Number:	206752-36-5, 1670-14-0
Chemical Formula:	$C_7H_8N_2 \cdot HCl \cdot H_2O$
Molecular Weight:	174.63

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Assay (HPLC):	min 99%
Assay (Titration):	min 98%
Water content (Karl Fischer):	max 12%
Solubility:	Clear to slightly hazy, colourless to faint yellow