



SPECIFICATION

**2-Acetamido-3,4,6-tri-O-acetyl-2-deoxy- α -D-glucopyranosyl
chloride
Stabilised with 2.5% CaCO₃**

Product Code:	MA06734
Synonyms:	2-Acetamido-3,4,6-tri-O-acetyl-2-deoxy- α -D-glucopyranosyl chloride 2-Acetamido-2-deoxy- α -D-glucopyranosyl chloride 3,4,6-triacetate
CAS Number:	3068-34-6
Chemical Formula:	C ₁₄ H ₂₀ ClNO ₈
Molecular Weight:	365.76

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 95%
Specific optical rotation [α] _D	+100 to +120° ²⁰ (c=1, CHCl ₃):