



## SPECIFICATION

### 5-Bromo-4-chloro-3-indolyl b-D-glucuronide cyclohexylammonium salt

Product Code:	EB04274
Synonyms:	X-Glucuronide X-Glc CHA-X-Glc
CAS Number:	114162-64-0
Chemical Formula:	$C_{14}H_{13}BrClNO_7 \cdot C_6H_{13}N$
Molecular Weight:	521.79

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Assay:	Min 99%
Purity (HPLC):	Min 99%
Identity (IR):	Conforms to structure
Specific Optical Rotation [ $\alpha$ ] <sub>D</sub>	- 87.0 to -91.0° <sup>25</sup> (c=1, DMF:H <sub>2</sub> O, 1:1):
Water content (K.Fischer):	2.5 - 4.5%