



## SPECIFICATION

### N6-Benzoyl-2'-deoxy-5'-O-DMT-adenosine 3'-CE phosphoramidite

Product Code:	PB09172
CAS Number:	98796-53-3
Chemical Formula:	$C_{47}H_{52}N_7O_7P$
Molecular Weight:	857.93

#### TECHNICAL SPECIFICATION

Appearance:	White to pale yellow crystalline powder
Purity (HPLC):	min 98%
Identity ( $^{31}P$ NMR):	Conforms to structure