



SPECIFICATION

6-Chloro-3-indolyl 2-acetamido-2-deoxy-b-D-glucopyranoside

Product Code:	EC05996
Synonyms:	6-Chloro-3-indolyl N-acetyl-b-D-glucosaminide 6-Chloro-3-(2-acetamido-2-deoxy-b-D-glucopyranosyloxy)indole
CAS Number:	156117-44-1
Chemical Formula:	$C_{16}H_{19}ClN_2O_6$
Molecular Weight:	370.78

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 97%