



## SPECIFICATION

### 6-Chloro-3-indolyl b-D-glucopyranoside

Product Code:	EC03297
Synonyms:	Rose glucoside
CAS Number:	159954-28-6
Chemical Formula:	$C_{14}H_{16}ClNO_6$
Molecular Weight:	329.73

#### TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 98%
Identity (IR):	Conforms to structure
Water content (Karl Fischer):	max 15%