



SPECIFICATION

2,2,6,6-Tetramethylpiperidine 1-oxyl, free radical - solid melt

Product Code:	FT04400
Synonyms:	TEMPO Tempol 4-Oxypiperidol
CAS Number:	2564-83-2
Chemical Formula:	$C_9H_{18}NO$
Molecular Weight:	156.25

TECHNICAL SPECIFICATION

Appearance:	Orange to red solid
Assay (GC):	min 98%
Identity (IR):	Conforms to structure
Residual solvents (as methanol):	max 0.3%