



## SPECIFICATION

### Lactodifucotetraose

|                   |                      |
|-------------------|----------------------|
| Product Code:     | OL06567              |
| Synonyms:         | LDFT                 |
| CAS Number:       | 20768-11-0           |
| Chemical Formula: | $C_{24}H_{42}O_{19}$ |
| Molecular Weight: | 634.58               |

## TECHNICAL SPECIFICATION

|                        |                       |
|------------------------|-----------------------|
| Appearance:            | White powder          |
| Purity:                | min 90%               |
| Identity ( $^1H$ NMR): | Conforms to structure |