



SPECIFICATION

1,3,4,6-Tetra-O-acetyl-2-azido-2-deoxy-b-D-glucopyranose

Product Code:	MT06641
CAS Number:	80321-89-7
Chemical Formula:	$C_{14}H_{19}N_3O_9$
Molecular Weight:	373.32

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 97%
Specific Optical Rotation $[\alpha]_D^{20}$	+2 to +4°

²⁰ (c=1, CHCl₃):