



SPECIFICATION

Sucrose monocaprate

Product Code:	DS06361
Synonyms:	b-D-Fructofuranosyl monodecanoate a-D-glucopyranoside n-Decanoylsucrose Sucrose monodecanoate
CAS Number:	31835-06-0
Chemical Formula:	$C_{22}H_{40}O_{12}$
Molecular Weight:	496.55

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 97%
Residue on Ignition:	max 1%
Water content:	max 1%
NOTE:	This product is a mixture of esters