



SPECIFICATION

5-Bromo-4-chloro-3-indolyl N-acetyl-a-D-N-acetylneuraminic acid sodium salt

Product Code:	EB02506
Synonyms:	X-NANA.Na
CAS Number:	160369-85-7
Chemical Formula:	$C_{19}H_{21}BrClN_2O_9.Na$
Molecular Weight:	559.72

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 90%
Solubility (50mg/ml in water):	Clear, colourless to pale yellow solution
Water content (Karl Fischer):	max 9.5%