



SPECIFICATION

Digoxin

Product Code:	OD01513
CAS Number:	20830-75-5
Chemical Formula:	$C_{41}H_{64}O_{14}$
Molecular Weight:	780.94

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity:	min 95%
Specific optical rotation:	+10 to +13°
Loss on drying:	max 1%
Identity (IR):	Conforms to structure
Solubility (50 mg/ml in pyridine):	Clear to very slightly hazy colourless to faint yellow solution
Sulfated ash:	max 0.1%