



SPECIFICATION

2-Acetamido-2-deoxy-6-O-(α -L-fucopyranosyl)-D-glucopyranose

Product Code:	OA04444
Synonyms:	Fuc-a-1,6-GlcNAc
CAS Number:	33639-80-4
Chemical Formula:	$C_{14}H_{25}NO_{10}$
Molecular Weight:	367.35

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure