



SPECIFICATION

5-Bromo-4-chloro-3-indolyl α -D-glucopyranoside

Product Code:	EB05081
Synonyms:	X-a-D-Glucoside 5-Bromo-4-chloro-3-(α -D-glucopyranosyloxy)indole
CAS Number:	108789-36-2
Chemical Formula:	$C_{14}H_{15}BrClNO_6$
Molecular Weight:	408.63

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 97%