



SPECIFICATION

5-Bromo-4-chloro-3-indolyl b-D-cellobioside

Product Code:	EB03309
CAS Number:	177966-52-8
Chemical Formula:	$C_{20}H_{25}BrClNO_{11}$
Molecular Weight:	570.77

TECHNICAL SPECIFICATION

Appearance:	White crystalline powder
Purity (HPLC):	min 99%
Specific optical rotation $[\alpha]_D^{20}$	-66.5 to -69.5° ²⁰ (c=1, MeOH):
Water content (Karl Fischer):	max 5.0%