

## PRODUCT SPECIFICATION



*Cat. No.*     **Z-4005\_P00**  
*Product*     **STAADIUM(TM) PeptiZide Pyr, Patent pending**

*CAS RN*

*M. F.*    C<sub>12</sub>H<sub>14</sub>N<sub>2</sub>O<sub>3</sub>

*M. W.*    409.41 g/mol

*Test*

*Specification*

**Color**

white to slightly red

**Form**

powder

**Purity (TLC)**

min. 95 area-%

(20 mg/mL in methanol/water 1:1; 1-butanol/acetic acid/water  
6:1:1; UV 254 nm, visual)

**<sup>1</sup>H-NMR Spectrum**

corresponds to assigned structure

(DMSO)

**Recommended Retest Period:**

**720 days**

*Date of issue:* October 6, 2017

**BIOSYNTH AG**

Luca Deriu