

## PRODUCT SPECIFICATION



*Cat. No.* **F-6660**

*Product* **6-Fluorotryptamine hydrochloride**

**3-(2-Aminoethyl)-6-fluoroindole hydrochloride**

*CAS RN* [55206-24-1]

*M. F.* C<sub>10</sub>H<sub>12</sub>ClFN<sub>2</sub>

*M. W.* 214.67 g/mol

<i>Test</i>	<i>Specification</i>
<b>Color</b>	white to off-white
<b>Form</b>	crystalline powder
<b>Water Content (K. Fischer)</b>	max. 1.0 %
<b>Quantitative 19F-NMR</b> (calculated on dry matter)	min. 98 %
<b>FT-IR Spectrum</b>	corresponds to reference spectrum
<b>1H-NMR Spectrum</b>	corresponds to assigned structure
<b>13C-NMR Spectrum</b>	corresponds to assigned structure
<b>19F-NMR Spectrum</b>	corresponds to assigned structure
<b>Recommended Retest Period:</b>	<b>1440 days</b>

*Date of issue:* March 6, 2019

**BIOSYNTH AG**

Susanne Babirat